



Experience of Hybrid Planning Technique in Breast Cancer Radiotherapy







Delta Hospital Limited

Md. Abu Kausar

Senior Medical Physicist

Dept. of Clinical Oncology & Medical Physics

Key Focuses

-  **Basics of Hybrid VMAT**
-  **Comparison of 3DCRT, VMAT and Hybrid Planning**
-  **Clinical benefits, dosimetric considerations and outcomes**
-  **Summary**

Breast Cancer Quick Facts

GLOBOCAN cancer incidence and mortality for 36 cancers in 185 countries





TABLE 1 New cases and deaths for 36 cancers and all cancers combined in 2022.

Cancer site	Incidence			Mortality		
	Rank	New cases	% of all sites	Rank	Deaths	% of all sites
Lung	1	2,480,301	12.4	1	1,817,172	18.7
Female breast	2	2,308,897	11.6	2	1,665,684	6.9
Colorectum	3	1,926,118	9.6	2	903,859	9.3
Prostate	4	1,466,680	7.3	8	396,792	4.1

- In 2024- more than 820 patients received radiotherapy including 2D, 3DCRT & IMRT/VMAT

Year	Right Breast Patients	Left Breast Patients
2024	392	428

Radiotherapy Left Breast Cancer

Multimodal Treatment for Locally Advanced Breast Cancer (LABC)



Standard Treatment Approach:

- Modified radical mastectomy followed by chemotherapy and radiotherapy.



Benefits of Post-Mastectomy Radiotherapy (PMRT):

- **M Clarker et al. Lancet 2005;366(9503):2087-106**
 - Reduces **chest wall (CW) recurrence** by a factor of **3–4** and improves **overall survival** by **6%**.

Radiotherapy Left Breast Cancer





Radiation Therapy Techniques in Breast Cancer:

- 3DCRT, IMRT, and VMAT are commonly used approaches.



3DCRT: The Reference Technique

-  Widely available and cost-effective.
-  Challenges:
- Low conformity & homogeneity.
- Under-coverage at the **CW-SCL junction**.
- High-dose spillage outside the **Planning Target Volume (PTV)**.


Radiotherapy Left Breast Cancer



- **IMRT & VMAT help**
 - to reduce dose to the heart and IL
 - improving the PTV coverage, conformity and homogeneity
 - VMAT has the advantage of a faster treatment time.
- **However, VMAT in breast cancer**
 - induces a considerable increase in the volume of healthy tissue receiving low doses
 - thereby increasing the risk of secondary cancers

Why reconsider traditional VMAT?

- **Important to account for low-dose exposure to healthy tissue**
- **Byrne M et al. J Appl Clin Med Phys 2018; 19: 684–93.**
 - **The generation of low-dose baths specifically to the contralateral breast, lung and heart which exceed that of 3DCRT**
- **Ma et al. Sci Rep 2015; 5: 12274**
 - **VMAT was associated with an increase in mean heart dose and low-dose volume to the lung, compared with 3D-CRT, possibly explaining the higher use of tangential IMRT over VMAT**
- **Haciislamoglu et al. BJR 2019; 92: 20190317**
 - **An estimation of excess absolute risk (EAR) in terms of developing a secondary cancer and observed the volume of low dose (3 Gy and 5 Gy) to normal tissue was significantly higher with IMRT and VMAT than FiF**



Solution – Hybrid VMAT



To improve outcomes, we implemented:

- **HYBRID VMAT** approach
- Combining VMAT with the **3D Field-in-Field (FiF) technique**
- Aiming to optimize target coverage while reducing unnecessary exposure to healthy tissue

Some Background on Hybrid VMAT...

www.nature.com/scientificreports

scientific reports

OPEN

Hybrid VMAT-3DCRT as breast cancer treatment improvement tool

Cyril Voyant^{1,2,3}, Morgane Pinpin², Delphine Leschi², Séverine Prapant²,
Françoise Savigny² & Marie-Aimée Acquaviva²


Radiation therapy is an important tool in the treatment of breast cancer and can play a crucial role in improving patient outcomes. For breast cancer, if the technique has been for a long time the use of 3DCRT, clinicians have seen the management evolve greatly in recent years. Field-in-field and IMRT approaches and more recently dynamic arthrotherapy are increasingly available. All of these approaches are constantly trying to improve tumour coverage and to preserve organs at risk by minimising the doses delivered to them. If arthrotherapy allows a considerable reduction of high doses received by healthy tissues, no one can deny that it also leads to an increase of low doses in tissues that would not have received any with other techniques. We propose a hybrid approach combining the robustness of the 3DCRT approach and the high technicality and efficiency of arthrotherapy. Statistical tests (ANOVA, Wilcoxon, determination coefficient, ROC, etc.) allow us to draw conclusions about the possibility of using the hybrid approach in certain cases (right breast, BMI > 23, age > 48, target volume > 350 cc, etc.). Depending on the breast laterality and patients morphological characteristics, hybridization may prove to be a therapeutic tool of choice in the management of breast cancer in radiotherapy.

Medical Dosimetry 40 (2015) 262–267



Medical Dosimetry

journal homepage: www.meddos.org



Dosimetric comparison of hybrid volumetric-modulated arc therapy, volumetric-modulated arc therapy, and intensity-modulated radiation therapy for left-sided early breast cancer

Jia-Fu Lin, M.S.^{*}, Dah-Cherng Yeh, M.D.[†], Hui-Ling Yeh, M.D.[‡], Chen-Fa Chang, M.S.^{*} and Jin-Ching Lin, M.D., Ph.D.[‡]

^{*}Department of Radiation Physics, Taichung Veterans General Hospital, Taichung, Taiwan; [†]Department of General Surgery, Taichung Veterans General Hospital, Taichung, Taiwan; and [‡]Department of Radiation Oncology, Taichung Veterans General Hospital, Taichung, Taiwan



Physica Medica 52 (2018) 86–92



ELSEVIER

Contents lists available at ScienceDirect

Physica Medica

journal homepage: www.elsevier.com/locate/ejmp



Original paper

Hybrid volumetric modulated arc therapy for chest wall irradiation: For a good plan, get the right mixture

Karunakaran Balaji^{a,b,*}, Poonam Yadav^c, Sitaraman BalajiSubramanian^a,
Chandrasekaran Anu Radha^b, Velayudham Ramasubramanian^b

^aDepartment of Radiation Oncology, Gleneagles Global Hospitals, Chennai, India

^bSchool of Advanced Sciences, Vellore Institute of Technology, Vellore, India

^cDepartment of Human Oncology, School of Medicine and Public Health, University of Wisconsin-Madison, Madison, WI, United States



Background...



ROJ Radiation
Oncology
Journal

Original Article

pISSN 2234-1900 · eISSN 2234-3156
Radiat Oncol J 2020;38(4):270-281
<https://doi.org/10.3857/roj.2020.00619>

Dosimetric comparative study of 3DCRT, IMRT, VMAT, Ecomp, and Hybrid techniques for breast radiation therapy

Semaya Natalia Chen^{1,2}, Prabhakar Ramachandran^{1,2}, Pradip Deb¹

¹School of Health and Biomedical Sciences, RMIT University, Melbourne, Australia

²Department of Radiation Oncology, National Cancer Centre Singapore, Singapore

³Princess Alexandra Hospital, Queensland, Australia

La radiologia medica
<https://doi.org/10.1007/s11547-019-00994-1>

MEDICAL PHYSICS



Hybrid volumetric modulated arc therapy for whole breast irradiation: a dosimetric comparison of different arc designs

Velayudham Ramasubramanian¹ · Karunakaran Balaji^{1,2} · Sitaraman Balaji Subramanian² · Krishnamoorthi Sathya² · Moorthi Thirunavukarasu² · Chandrasekaran Anu Radha¹

Received: 8 June 2018 / Accepted: 24 January 2019
© Italian Society of Medical Radiology 2019

Journal of Radiation Research, Vol. 61, No. 5, 2020, pp. 747–754
doi: 10.1093/jrr/rraa057
Advance Access Publication: 12 August 2020



Hybrid volumetric-modulated arc therapy for postoperative breast cancer including regional lymph nodes: the advantage of dosimetric data and safety of toxicities

Yoshiko Doi^{1,2,*}, Minoru Nakao¹, Hideharu Miura¹, Shuichi Ozawa¹,
Masahiro Kenjo¹ and Yasushi Nagata^{1,2}

¹Hiroshima High-Precision Radiotherapy Cancer Center, 3-2-2, Futatabanosato, Higashi-ku, Hiroshima-shi, Hiroshima, 732-0057, Japan

²Department of Radiation Oncology, Institute of Biomedical & Health Sciences, Hiroshima University, 3-2-2, Futatabanosato, Higashi-ku, Hiroshima-shi, Hiroshima, 732-0057, Japan

*Corresponding author. Hiroshima High-Precision Radiotherapy Cancer Center, 3-2-2, Futatabanosato, Higashi-ku, Hiroshima-shi, Hiroshima, 732-0057, Japan.

Tel: +81-82-263-1330; Fax: +81-82-263-1331; Email: doi@hiprac.jp

(Received 17 February 2020; revised 27 May 2020; editorial decision 29 June 2020)


Workflow and Implementation

❖ Key Considerations:

- 10 Lt Breast Patients
- Chest wall + Suprclav
- Fraction size/dose
 - * 2.68 Gy/15 fractions

☐ We delineate

- Heart
- Ipsilateral Lung
- Contralateral Lung
- Contralateral Breast
- Chestwall PTV volume



Workflow and Implementation

❖ Planning consideration:

- Varian Eclipse Treatment Planning System
- True Beam Linear Accelerator
- 3DCRT
- VMAT
- Hybrid (3DCRT+VMAT)

For Each Patients we established

- 6 dosimetrics
- 3DCRT
- VMAT

Hybrid

80% 3DCRT/20% VMAT

70% 3DCRT/30% VMAT

60% 3DCRT/40% VMAT


50% 3DCRT/ 50% VMAT

Treatment Planning

- Isocentre positioned between CW and SCL junction
- Same isocentre for VMAT, 3DCRT and hybrid
- Same Collimator angle (0°)
- Gantry angle selected based on PTV curvature

Method	Prescription	Beams	Nature	Gantry Angle	Rotation	Colli
3DCRT	CW+SCL	Med. Tangent	FiF/6X	305-310	n/a	0
		Lat. Tangent	FiF/6X	125-130	n/a	0
VMAT	CW+SCL	Arc 1	VMAT/6X	305-310	CW	0
		Arc 2	VMAT/6X	135-140	CCW	0
Hybrid	CW+SCL	Tangent+Arc	FiF+VMAT /6X	3D+VMAT	n/a & Arc	0

Treatment Planning.....




Planning Objectives

PTV/OARs	Dose Constraints
PTV	$D95\% \geq 38.19 \text{ Gy}$ $D2\% \leq 44.22 \text{ Gy}$ $HI = 0$ $CI \geq 0.80$
Heart	$V20\text{Gy} \leq 5 \%$ $\text{Mean} \leq 4 \text{ Gy}$
Ipsilateral Lung	$V20\text{Gy} \leq 15 \%$ $V10\text{Gy} \leq 35 \%$ $V5\text{Gy} \leq 50 \%$ $\text{Mean} \leq 12 \text{ Gy}$
Contralateral Lung	$\text{Mean} \leq 3\text{Gy}$
Contralateral Breast	$\text{Mean} \leq 3\text{Gy}$

RADIATION THERAPY ONCOLOGY GROUP

RTOG 1005


A PHASE III TRIAL OF ACCELERATED WHOLE BREAST IRRADIATION
WITH HYPOFRACTIONATION PLUS CONCURRENT BOOST
VERSUS STANDARD WHOLE BREAST IRRADIATION
PLUS SEQUENTIAL BOOST FOR EARLY-STAGE BREAST CANCER





**KEEP
CALM
AND
BEAM
ON**

Comparison of 3DCRT & VMAT


Parameter		3DCRT	VMAT
PTV	D95	37.70	39.00
	D2	43.50	43.20
	HI	0.20	0.13
	CI	1.50	1.04
Heart	Mean	3.50	5.10
	V20	6.00	2.00
Lung Ipsi	V20	23.00	21.00
	V10	30.10	40.30
	V5	42.10	69.5
	Mean	11.0	13.1
Lung Contra	Mean	0.20	3.40
Breast Contra	Mean	0.30	3.60



3CRT




VMAT




Different Proportion of Hybrid VMAT

Parameter		80%/20% H-VMAT	70%/30% H-VMAT	60%/40% H-VMAT	50%/50% H-VMAT
PTV	D95	38.50	38.60	38.60	38.70
	D2	43.20	43.0	42.9	43.10
	HI	0.15	0.13	0.13	0.14
	CI	1.37	1.20	1.15	1.13
Heart	Mean	3.80	3.83	4.10	4.53
	V20	5.30	5.31	4.96	4.65
Lung Ipsi	V20	23.60	23.2	22.8	22.4
	V10	31.8	32.40	33.30	34.6
	V5	47.4	49.8	52.30	55.8
	Mean	11.30	11.15	11.38	11.78
Lung Contra	Mean	0.8	1.20	1.50	1.80
Breast Contra	Mean	0.95	1.25	1.60	1.90




3DCRT, VMAT & Hybrid





3DCRT & Hybrid



VMAT & Hybrid



Comparative dose distribution of different dose weighting




Summary of Clinical Studies on Whole-Breast Radiotherapy (WBRT) Using Hybrid Techniques.

Study	Dose Weighting	Findings
Abrose et al.(2021)	60%–80% 3D + 20%–40% IMRT	H-IMRT provided a good dose coverage and OAR protection.
Cilla et al.(2021)	80% 3D-CRT + 20% VMAT	H-VMAT resulted in superior target dose conformity and homogeneity compared to other techniques.
Cunningham et al.(2021)	70% 3DCRT + 30% IMRT	H-IMRT and IMPT provided better homogeneous dose distribution than IMRT.
Chen et al. (2020)	70% 3D CRT + 30% VMAT	H-VMAT advantages in the conformity index. H-VMAT plans better to the H-IMRT plans with regard to heart dose and treatment delivery time.
Balaji et al (2020)	70% 3D-CRT + 30% VMAT	H-IMRT and H-VMAT provide similar dosimetric results.
Chan et al.(2017)	80% 3D-CRT 20% IMRT	H-VMAT advantages in the conformity index. H-VMAT plans better to the H-IMRT plans with regard to heart dose and treatment delivery time.

Do we need Skin Flash?

- Open fields typically contribute 70-80% of the total weight and include flash to compensate for positioning, breathing motion, and setup errors.
- Maintaining this ratio may eliminate the need for flash, but reducing it requires retaining flash.



Summary



- A **70–80% dose weighting** from **3DCRT** and **20–30% from VMAT** was found to be optimal for the **H-VMAT technique**.
- If hybrid planning is the preferred approach for target volume coverage, **3DCRT for OAR sparing**.
- Lung CL and LungIL, there is an extremely **significant difference highlighting** the quality of the 3DCRT & VMAT.
- In this study, we realize that the doses to the **contralateral breast and lung are indeed low ($D_{\text{mean}} < 3 \text{ Gy}$)** but much higher than those observed with 3DCRT.

- It is important to mention that the increase of **low dose bath** will provide in principle an **increase in second cancer risk**, at different years post radiotherapy.

Thank You

IF I CAN DO
CANCER
I CAN DO
ANYTHING

